

# DermaPure Wound and Skin Cleanser

SNS51118-8oz. bottle

## Instructions for Use

# **DermaPure Antimicrobial Wound and Skin Cleanser**

Simpurity<sup>™</sup> DermaPure Antimicrobial Wound and Skin Cleanser is recommended by doctors for the care of chronic wounds. DermaPure wound cleanser is steroid and alcohol free and contains no fragrance, dyes or perfumes. The formula is hypoallergenic and non-irritating to the skin and wound bed. Contains .018% Hypochlorous Acid.

### **Application:**

1. Spray the DermaPure cleanser onto the wound bed to cleanse the wound 2. Wipe with a sterile gauze cleaning sponge

### Indications:

Simpurity<sup>™</sup> DermaPure Wound and Skin Cleanser is indicated for the management of wounds such as: minor burns, skin abrasion, superficial wounds, minor lacerations, palliative wounds, bed sores, sloughy, necrotic and infected wounds.

#### **Contraindications:**

The product is not indicated for patients with a known sensitivity to hypochlorous acid and the cleanser's other ingredients.