



Product Overview

Safe n Simple focuses on product quality, patient safety, and cost containment. With the introduction of the Simpurity Silver Foam wound care dressing we are proud to expand the Simpurity wound care product line. Simpurity wound care products follow evidence based medicine, facilitate moist wound healing and improve patient outcomes while lowering costs.

Simpurity Silver Foam Dressing is a silver ion-containing polyurethane foam dressing that protects wounds from bacterial penetration while reducing the bioburden in a wound bed. It mainly consists of a hydrophilic polyurethane foam substrate and silver compound. It is indicated for the management of moderate-to-heavy exudate for chronic and acute wounds and helps maintain a moist healing environment. Controlled release of ionic silver provides a sustained antibacterial effect. Inhibits a broad-spectrum of infection up to 99.99% for a wide variety of bacteria. Provides longer wear times and pain-free removal.

Simpurity Silver Foam Dressings are available in 2"x2" and 4"x5" sizes.





Product Features

- Sterile, foam dressing helps to maintain an optimal balance in fluid creating a moist wound environment for healing
- Easy to apply and remove
- High absorption capacity
- Protective backing to protect against strikethrough
- Impregnated with ionic silver particles to kill bacteria

Indications

Simpurity Silver Foam Wound Dressings are indicated for the management of moderately-to-highly exuding wounds such as:

- Pressure ulcers
- Leg ulcers
- Diabetic foot ulcer
- Second degree burns
- Donor sites
- Post operative wounds and skin abrasions

Product Details

Product	Item Number	Units/Box	Boxes/Case	HCPC
2" x 2"	SNS72502	12	16	A6209
4" x 5"	SNS72520	12	16	A6210
4" x 6"	SNS72524	12	16	A6210
6" x 6"	SNS72506	12	16	A6210
8" x 8"	SNS72508	12	16	A6211

