



## Product Overview

Safe n Simple focuses on product quality, patient safety, and cost containment. With the introduction of the Simpurity Hydrocolloid wound care dressing we are proud to expand the Simpurity wound care product line. Simpurity wound care products follow evidence based medicine, facilitate moist wound healing and improve patient outcomes while lowering costs.

Simpurity Hydrocolloid wound dressings accelerates the healing process by providing an ideal environment for formation of granulation tissue. The dressing forms a gel upon absorbing wound exudates providing a moist wound environment, supporting granulation and epithelialization. The gel formation prevents Simpurity Hydrocolloid from adhering to wounds and dressing can be removed without causing pain and without damaging the granulation and epithelial tissue.

Simpurity Hydrocolloid is a gelatin free hydrocolloid adhesive wound dressing laminated to a semi permeable, bacteria proof, virus proof and waterproof polyurethane film backing individually packed and sterilized by radiation. Simpurity Hydrocolloids high absorption capacity binds away excess exudate prolonging dressing change intervals. The dressing is slightly transparent enabling wound inspection.



## Product Features

- Delivers gentle and secure adhesion that does not irritate the skin
- Provides a barrier to external debris and bacteria when securely adhered to the skin
- Has a silicone release liner to facilitate removal from adhesive and application to the skin
- Latex free

## Indications

Simpurity Hydrocolloid Wound Dressings are indicated for the management of light-to-moderately exuding wounds such as:

- Second-degree burns
- Leg ulcers
- Pressure ulcers
- Chronic slow healing wounds and abrasions
- Used also for blister treatments

## Contra-Indications

Simpurity Hydrocolloid Wound Dressings are not indicated for use on the following

- Third-degree burns
- On exposed muscle, bones and tendons
- Infected wounds

## Product Details - Hydrocolloid

Product	Item Number	Units/Box	Boxes/Case	HCPC
2" x 2"	SNS55402	12	12	A6234
4" x 4"	SNS55404	12	12	A6235



Clarkston, MI 48348 | 844-767-6334 | [www.SnS-Medical.com](http://www.SnS-Medical.com)

WCF-030116