Introducing

Silicone

Foam with Adhesive Border

Safe n Simple focuses on product quality, patient safety, and cost containment. Simpurity wound care dressings follow evidence based medicine, facilitates moist wound healing and improves patients outcomes while lowering costs.



Product Details:

Size	ltem#	Box/Case	HCPC
3"x3"	SNS77733	16 bxs/cs	A6212
4"x4"	SNS77744	16 bxs/cs	A6212

FEATURES:

- Core silicone oil technology provides super gentle silicone contact layer for fragile skin
- Highly absorbent foam with a gentle adhesive silicone contact layer and adhesive border to secure the dressing in place
- Latex Free
- Maintains a moist wound healing environment
- No irritation to the skin upon removal of dressing
- · Minimizes trauma to the wound

INDICATIONS:

For the management of moderate-to-highly exuding wounds such as:

- Pressure Ulcers
- Leg Ulcers
- · Non-infected diabetic foot ulcer
- Donor sites
- Post-operative wound
- Lacerations and abrasions

CLINICAL DATA:

- Contains a PU foam pad that can absorb more than 5 times its original pad weight
- Contains a silicone gel adhesive that is applied to be greater than 100g/sq meter
- Contains super absorbent fibers that is greater than 60 grams /sq meter
- Contains a polyurethane film to prevent strikethrough
- Contains a nonwoven fabric
- Contains a PE release film



- 180 degree peel strength is greater than .3N/25cm @ room temperature
- Fluid handling capacity is greater than 5000g/sq meter/24hrs
 @ 37 degrees C
- Absorption Greater than 5 times the dressing weight

