



# Wound Care Product Guide

# Introduction...

Safe n Simple™ is a medical care innovator, specializing in the promotion of healthy skin, complete security and confident ostomy and wound care management. We view our product development through the eyes of our customers, who want active, productive lives unimpeded by their conditions. Our goal is to make our customer's health care needs safe and simple.

We follow two overriding principles: create innovative, high-performance products, and deliver the greatest possible value. We continually push our technicians to develop better ideas, more advanced formulations and superior products, while at the same time making those products convenient and affordable.

And now with Simpurity Wound Dressings we offer a complete brand of wound care products created to promote healing and cost savings to our customers.

Above all, we provide our customers with safety, security and simplicity. With Safe n Simple, that's our promise.

# Table of Contents

Page 2-4 .....	Wound Care Intro Pages & Table of Contents
Page 5 .....	Alginate
.....	Silver Alginate
Page 6 .....	Collagen
.....	Fibergel
Page 7 .....	Silver Fibergel
.....	Borderless Foam
Page 8 .....	Bordered Foam
.....	Bordered Fenestrated Circular Foam
Page 9 .....	Ionic Silver Foam
.....	Borderless Silver Silicone Foam
Page 10.....	Borderless Silicone Foam
.....	Bordered Silicone Sacrum Foam
Page 11.....	Bordered Silver Silicone Foam
.....	Bordered Silicone Foam
Page 12.....	Borderless Gauze Sponges
.....	Bordered Gauze Sponges
Page 13.....	Microroll & Microroll Slim Gauzes
.....	Hydrocolloids
Page 14.....	Hydrogel Sheets
.....	Bordered Hydrogel Sheets
Page 15.....	Hydrogel Impregnated Gauze
.....	Silicone Transparent Film
Page 16.....	IV Derm Silicone Transparent Film
.....	Polymer Pads
Page 17.....	Border & Borderless Super Absorbent Pads
.....	DermaPro Waterproof Plastic Tape
Page 18.....	DermaPro Waterproof Silicone Tape
.....	Xeroform Impregnated Gauze
Page 19.....	DermaPure Antimicrobial Wound and Skin Cleanser
.....	Wound Care Allowables
Page 20 & 21.....	Wound Care Kits
Page 22-23 .....	Stay Connected with SNS



## Enhancing quality of life with exceptional products at an affordable price

Safe n Simple focuses on product quality, patient safety, and cost containment. With the introduction of the Simpurity Wound Care product line we are proud to offer wound care products that follow evidence based medicine, facilitate moist wound healing and improve patient outcomes while lowering costs.

Your patients deserve the best care available, with the proven quality of Simpurity wound care product solutions, you can help them achieve greater outcomes and save money doing it.

### Testimonial.....

"I personally use Safe n Simple wound care products in our wound center for several reasons. The Simpurity line of products offers high quality alginate, foam and collagen dressings with and without silver at a very cost effective price. When everyone is looking to use evidenced based treatments and work within our current reimbursement environment this is the PERFECT product line. The products are user friendly. The packaging is appealing and instructions clear. Most importantly our healing rates and results have been excellent."

**J. Christopher Brandys, MD FRCS FACS**  
*General and Vascular Surgeon  
Director of The Wound Center  
Medical Consultant - Safe n Simple*

For more information or to request samples please call our toll-free number or visit our website.

**844-767-6334** | [www.SnS-Medical.com](http://www.SnS-Medical.com)

Affordable products covered by insurance.

# Alginate

## Alginate Wound Dressing

Simpurity Alginate Dressings are derived from seaweed. The products are composed of soft, calcium alginate fibers shaped as pads (fibrous mats) or ropes (twisted fibers). When packed in a wound, the alginate interacts with wound exudate to form a soft gel that maintains a moist healing environment. The dressings release of calcium ions can activate prothrombin and improve hemostasis. Has a strong tensile strength.




Item	Product#	Size	HCPC
Dressing Pad, 10/Box	SNS50702	2"x2"	A6196
Dressing Pad, 10/Box	SNS50704	4"x4"	A6196
Dressing Pad, 10/Box	SNS50720	4"x5"	A6197
Dressing Pad, 5/Box	SNS50732	4"x8"	A6197
Dressing Pad, 10/Box	SNS50764	8"x8"	A6198
Rope, 30/Bag	SNS50712	12"	A6199



## Silver Alginate Wound Dressing

Simpurity Silver Alginate Dressings are made of natural fibers derived from seaweed through a series of special processes. Simpurity Silver Alginate Dressings can inhibit bacteria growing in the wound bed reducing risk of infection while maintaining a moist environment for wound healing. The dressing has a strong tensile strength, a high absorption capacity and is easy to remove from the wound bed.



Item	Product#	Size	HCPC
Dressing Pad, 10/Box	SNS51702	2"x2"	A6196
Dressing Pad, 10/Box	SNS51704	4"x4"	A6196
Dressing Pad, 10/Box	SNS51720	4"x5"	A6197
Dressing Pad, 10/Box	SNS51764	8"x8"	A6198
Rope, 10/Box	SNS51712	12"	A6199



# Collagens

## Collagen Wound Dressing

Simpurity Collagen provide ingredients for tissue regeneration in phase 3 wounds and the 3D scaffold structure of Simpurity Collagen allows fibroblasts to migrate across the wound bed facilitating granulation, angiogenesis and tissue repair. Simpurity Collagen dressings are 100% pure Type I, non-bleached, native undigested bovine collagen. The dressings provide high exudate absorption and retention providing a moist wound healing environment. Easy application and superior conformity to the wound.



Item	Product#	Size	HCPC
Dressing Pad, 10/Box	SNS50002	2"x2"	A6021
Dressing Pad, 5/Box	SNS52222	2"x2"	A6021
Dressing Pad, 5/Box	SNS52244	4"x4"	A6021
Vial Particles, 5/Box	SNS5001G	1 Gram	A6010
Packet Particles, 10/Box	SNS5221G	1 Gram	A6010

Note: SNS50002 is a premium brand collagen wound dressing.



# Fibergels

## Fibergel Wound Dressing

Simpurity Fibergel Wound Dressing is a soft, absorbent material that forms into a gel when in contact with exudate. The gelling action creates an optimal environment for moist wound healing. It conforms easily to various wound shapes, supports debridement and can easily be removed in one piece with minimal damage to newly forming granulating tissues.



Item	Product#	Size	HCPC
Dressing Pad, 10/Box	SNS56702	2"x2"	A6196
Dressing Pad, 10/Box	SNS56704	4"x4"	A6196
Dressing Pad, 10/Box	SNS56745	4"x5"	A6197



# Fibergels

## Silver Fibergel Wound Dressing

Simpurity Fibergel Ag Wound Dressing has all the benefits of a hydrofiber dressing plus ionic silver to inhibit bacteria from growing in the wound bed, reducing the risk of infection. The gelling action of Fibergel AG creates an optimal environment for moist wound healing. It conforms easily to various wound shapes, supports debridement and can easily be removed in one piece with minimal damage to newly forming granulating tissues.



Item	Product#	Size	HCPC
Dressing Pad, 10/Box	SNS56712	2"x2"	A6196
Dressing Pad, 10/Box	SNS51744	4"x4"	A6196
Dressing Pad, 10/Box	SNS56716	4"x4.75"	A6197



# Foams

## Foam Non-Adhesive Wound Dressing

Simpurity Foam Wound Dressings are designed to help manage moderate-to-highly draining wounds. The dressing absorbs wound exudate, promotes moist wound healing, protects the wound from contaminants and permits moisture vapor exchange. These foam dressings have a protective colored film backing to protect against strike through.



Item	Product#	Pad Size	HCPC
Dressing Pad, 12/Box	SNS51W02	2"x2"	A6209
Dressing Pad, 12/Box	SNS51W03	3"x3"	A6209
Dressing Pad, 12/Box	SNS51W04	4"x4"	A6209
Dressing Pad, 12/Box	SNS51W24	4"x6"	A6210
Dressing Pad, 12/Box	SNS51W08	8"x8"	A6211



# Foams

## Foam Wound Dressing - Adhesive Border

Simpurity Foam wound dressings are composed of a soft and highly absorbent polyurethane foam pad to absorb wound drainage and a transparent adhesive border to hold the dressing in place. This island foam dressing helps maintain a moist environment, which accelerates wound healing. It is especially suitable for moderate to highly exuding wounds and has a protective backing to protect against strike through. The dressings are comfortable, easy to apply, and easy to remove from the wound.



Item	Product#	Pad Size	HCPC
Dressing Pad, 12/Box	SNS72302	2"x2"	A6212
Dressing Pad, 12/Box	SNS72320	4"x5"	A6213

Overall dressing sizes including borders: 4"x4" and 6"x7"



## Foam Wound Dressing - Fenestrated Round Border

Simpurity Foam Fenestrated Round Bordered dressing has a non-adhesive foam pad that absorbs large amounts of fluid. The dressing is fenestrated to keep tubing in place and helps to wick away moisture at the skin surface. The adhesive border secures the dressing in place to the skin surface.



Item	Product#	Pad Size	HCPC
Dressing Pad, 12/Box	SNS51C04	2" Circle	A6212

Overall dressing size including borders: 4" Round.





# Foams

## Foam Wound Dressing - Ionic Silver

Simpurity Foam Dressing is a silver ion-containing polyurethane foam dressing that protects wounds from bacterial penetration while reducing the bioburden in a wound bed. It mainly consists of hydrophilic polyurethane foam substrate and silver compound. It is indicated for the management of moderate-to-high exudate for chronic and acute wounds and helps maintain a moist healing environment. Controlled release of ionic silver provides sustained antibacterial effect. Inhibits a broad-spectrum of infection up to 99.99% for a wide variety of bacteria. Provides longer wear times and pain-free removal.



---

Item	Product#	Pad Size	HCPC
Dressing Pad, 12/Box	SNS72502	2"x2"	A6209
Dressing Pad, 12/Box	SNS72520	4"x5"	A6210

---



## Foam Wound Dressing - Silver Silicone

Simpurity Silicone Foam Dressing is a polyurethane foam pad containing silver ions that has a gentle adhesive silicone contact layer to protect wounds from bacterial penetration while reducing the bioburden in the wound bed. It is indicated for the management of moderate-to-high exudate for chronic and acute wounds and helps maintain a moist healing environment. Controlled release of ionic silver provides sustained antibacterial effect. Provides longer wear times and pain-free removal.



---

Item	Product#	Pad Size	HCPC
Dressing Pad, 12/Box	SNS74422	2"x2"	A6209
Dressing Pad, 12/Box	SNS74445	4"x5"	A6210
Dressing Pad, 12/Box	SNS74466	6"x6"	A6210
Dressing Pad, 12/Box	SNS74477	7"x7"	A6211

---



# Foams

## Foam Wound Dressing - Silicone

Simpurity Silicone Foam Dressing is a superior absorbent dressing for moderate to high wound exudate management. It contains a highly absorbent foam pad and gentle adhesive silicone contact layer. These dressings help maintain a moist wound healing environment with no irritation upon removal, minimizes trauma to the wound and reduces pain during dressing changes.



Item	Product#	Pad Size	HCPC
Dressing Pad, 12/Box	SNS74502	2"x2"	A6209
Dressing Pad, 12/Box	SNS74520	4"x5"	A6210



## Silicone Sacrum Foam Wound Dressing

Simpurity Sacrum Foam wound dressings are composed of a soft absorbent pad that enhances patient comfort and absorbs excess fluids. The gentle silicone border adheres to the surrounding skin, not the wound bed. Helps to minimize pain and trauma at dressing changes. Minimizes pressure ulcers of the sacral region.



Item	Product#	Pad Size	HCPC
Dressing Pad, 12/Box	SNS77777	5"x5"	A6213

Overall dressing sizes including border: 7"x7"



# Foams

## Foam Wound Dressing - Bordered Silver Silicone

Simpurity Silver Silicone Foam Bordered Dressing is a superior absorbent dressing containing silver ions for moderate to high wound exudate management. It contains a highly absorbent foam pad impregnated with silver ions and gentle adhesive silicone contact layer with an adhesive border to hold the dressing in place. These dressings help maintain a moist wound healing environment with no irritation upon removal, minimizes trauma to the wound and reduces pain during dressing changes.



Item	Product#	Pad Size	HCPC
Dressing Pad, 12/Box	SNS77111	1"x1"	A6212
Dressing Pad, 12/Box	SNS77133	3"x3"	A6212

Overall dressing sizes including borders: 2"x2" and 4.6"x4.6"



## Foam Wound Dressing - Bordered Silicone

Simpurity Silicone Foam Bordered Dressing is a superior absorbent dressing for moderate to high wound exudate management. It contains a highly absorbent foam pad and gentle adhesive silicone contact layer and an adhesive border to hold the dressing in place. These dressings help maintain a moist wound healing environment with no irritation upon removal, minimizes trauma to the wound and reduces pain during dressing changes.



Item	Product#	Pad Size	HCPC
Dressing Pad, 12/Box	SNS77733	3"x3"	A6212
Dressing Pad, 12/Box	SNS77744	4"x4"	A6212

Overall dressing sizes including borders: 5"x5" and 6"x6"



# Gauzes

## Gauze Sterile Sponge

Simpurity Gauze Sterile Sponges are ideal for general wound cleansing, debriding, minor prepping, covering and packing wounds. It is a woven sponge made of 12 ply 100% USP cotton. Contains 2 sponges per packet.



Item	Product#	Pad Size	HCPC
Dressing Pad, 2 Pcs/pkt, 25/Box	SNS54444	4"x4"	A6402



## Bordered Gauze Sterile Sponge

Simpurity Bordered Gauze Sterile Sponges are a non-adherent gauze sponge island dressing with an adhesive border. These woven sponges are made of 12ply 100% USP cotton for use as a primary or secondary dressing. The sponge absorbs exudate for a cleaner wound bed and the adhesive border protects the wound from outside contaminants.



Item	Product#	Pad Size	HCPC
Dressing Pad, 30/Box	SNS52316	2"x2"	A6219
Dressing Pads, 30/Box	SNS52324	4"x4.75"	A6220

Overall dressing sizes including borders: 4"x4" and 6"x6.75"



# Gauzes

## Microroll Gauze Wound Dressings

Simpurity Microroll Gauze is made of prewashed, 100% woven gauze that has been fluffed to provide for a more bulky cushion to protect wound areas. The 6 ply construction of the dressing provides for excellent absorbency, superior breathability and fast-wicking action. The Microroll gauze rolls are infused with PHMB which provide a broad range of antiseptic effectiveness against gram-positive and gram-negative microorganisms, including several multi-drug resistant strains including MRSA. Microroll gauze provides a lower cost method of treatment.

Microroll Gauze Slim provides a less bulky protection the the wound area with the same PHMB protection and 6ply construction.



Item	Product#	Size	HCPC
Dressing Roll, 100/Case	SNS53344	4.5"x4.1yds	A6266
Microroll Gauze Slim Dressing Roll, 100/Case	SNS53311	4.5"x4.1yds	A6266


# Hydrocolloids

## Hydrocolloid Wound Dressing

Simpurity Hydrocolloid Wound Dressing is a thin hydrocolloid sheet. Simpurity Hydrocolloids are absorbent, and delivers a gentle and secure adhesion that provides a barrier to external debris and bacteria when securely adhered to the skin. For use on light-to-moderately exuding wounds.



Item	Product#	Size	HCPC
Dressing Pad, 12/Box	SNS55402	2"x2"	A6234
Dressing Pad, 12/Box	SNS55404	4"x4"	A6234
Dressing Pad, 10/Box	SNS21149	7"x7"	A6236



# Hydrogels

## Hydrogel Dressing Sheet

Simpurity Hydrogel Absorbent Dressing Sheets are designed to provide a moist wound healing environment for wounds with minimal to no exudate. It is permeable to air and moisture vapor for faster wound healing. Easy to remove from the skin without damaging tissues. Conforms to the wound bed preventing bacterial infection and contamination.



Item	Product#	Pad Size	HCPC
Dressing Pad, 10/Box	SNS58302	2"x2"	A6242
Dressing Pad, 10/Box	SNS58345	4"x5"	A6243



## Hydrogel Bordered Dressing Sheet

Simpurity Hydrogel Absorbent Dressing Sheets have a polyurethane border and are designed to provide a moist wound healing environment for wounds with minimal to no exudate. It is permeable to air and moisture vapor for faster wound healing. Easy to remove from the skin without damaging tissues. Conforms to the wound bed preventing bacterial infection and contamination.



Item	Product#	Pad Size	HCPC
Dressing Pad, 12/Box	SNS58312	2"x2"	A6245
Dressing Pad, 12/Box	SNS58315	4"x5"	A6246

Overall dressing sizes including borders: 4"x4" and 6"x7"




# Hydrogels

## Hydrogel - Impregnated Gauze

Simpurity Hydrogel Impregnated Gauze Wound Dressing is designed to provide a moist wound healing environment for wounds with minimal or no exudate. Contains a high viscosity gel impregnated into the gauze pad. Easy to apply.



Item	Product#	Pad Size	HCPC
Dressing Pad, 10/Box	SNS58802	2"x2"	A6231
Dressing Pad, 10/Box	SNS58820	4"x5"	A6232



# Silicone Transparent Films

## Transparent Silicone Film Wound Dressing

Simpurity Transparent Silicone Film wound dressings provide a waterproof barrier to external contaminants. It is permeable to air and moisture vapor for faster wound healing. Gentle to remove from sensitive skin areas, can be changed without pain or damage to the skin.



Item	Product#	Size	HCPC
Dressing Sheet, 12/Box	SNS57245	4"x5"	A6258



# Silicone Transparent Films

## IV Derm Silicone Film Dressing

Simpurity IV Derm Transparent Silicone Film dressings provide a waterproof barrier to secure IV lines from outside contaminants. Gentle to remove from sensitive skin areas, can be changed without pain or damage to the skin.



---

Item	Product#	Size	HCPC
Dressing Sheet, 12/Box	SNS52723	2"x3"	A6257

---



# Super Absorbent Pads

## Polymer Wound Dressing - Highly Absorbent

Simpurity Polymer Foam Wound Dressing is a thin, sterile, superabsorbent single-use dressing, for use on highly exuding wounds. Simpurity Polymer Foam absorbs the exudate and locks it into the super-absorbent core, allowing for significantly more fluid to be handled by the dressing while keeping it away from the wound and surrounding tissue. This can help to reduce maceration, as well as reduce exposure to harmful components of wound fluid.



---

Item	Product#	Pad Size	HCPC
Dressing Pad, 10/Box	SNS59020	4"x5"	A6251
Dressing Pad, 10/Box	SNS59080	8"x10"	A6253

---





# Super Absorbent Pads

## Super Absorbent Pads Wound Dressing

Simpurity Super Absorbent Pads are made of a non-woven super absorbent fiber that will absorb up to 20x its weight of fluid. One layer of Super Absorbent Pad is equal to 36 layers of medical gauze. The pads will maintain a moist wound environment by allowing for breathability of moisture and oxygen for wound healing. Anti-seepage to prevent wound from leaking. A great alternative to ABD pads.



### Bordered Super Absorbent Pads

Item	Product#	Pad Size	HCPC
Dressing Pad, 10/Box	SNS59122	2"x2"	A6251
Dressing Pad, 10/Box	SNS59144	4"x4"	A6251

### Non-Bordered Super Absorbent Pads

Item	Product#	Pad Size	HCPC
Dressing Pad, 10/Box	SNS59022	2"x2"	A6251
Dressing Pad, 10/Box	SNS59044	4"x4"	A6251

Overall dressing sizes including borders: 4"x4" and 6"x6"



## Tapes

### DermaPro Waterproof Plastic Tape

Simpurity DermaPro Waterproof Plastic Tape is superior for securing wound dressings and repelling moisture from outside contaminants. DermaPro has a latex free adhesive and is gentle to remove from sensitive skin areas. The roll contains a start strip for easy dispensing.



Item	Product#	Size	HCPC
Plastic Tape, 1 Roll/Box	SNS57220	2"x15'	A4452

# Tapes

## DermaPro Waterproof Silicone Tape

Simpurity DermaPro Waterproof Silicone Tape is latex free and has a silicone adhesive that is gentle to remove from sensitive skin areas. Superior for securing wound dressings and repelling moisture from outside contaminants. Each roll is perforated for easy tearing. Contains 12-1" rolls per box or 6-2" rolls per box.



Item	Product#	Size	HCPC
Silicone Tape, 12 Rolls/Box	SNS57230	1"x15'	A4452
Silicone Tape, 6 Rolls/Box	SNS57232	2"x15'	A4452


### Want to try a sample of the silicone tape?

We have 2.5' rolls available in 1" and 2" widths  
Call our Customer Service  
to request yours today!

# Xeroforms

## Xeroform - Impregnated Gauze

Simpurity Xeroform Impregnated Gauze is a sterile, occlusive, non-adhering dressing consisting of an absorbent, fine-mesh gauze impregnated with an antimicrobial petrolatum blend. The dressing will conform to all body contours which makes it suitable for use in compression therapy, deep and shallow wounds, burns, plastic surgery and surgical wounds with light exudate where mild medication and deodorization are desired.



Item	Product#	Pad Size	HCPC
Dressing Pad, 10/Box	SNS58818	1"x8"	A6222
Dressing Pad, 10/Box	SNS58822	2"x2"	A6222
Dressing Pad, 10/Box	SNS58844	4"x4"	A6222
Dressing Pad, 10/Box	SNS58859	5"x9"	A6223



# Wound and Skin Cleanser

## DermaPure Antimicrobial Wound and Skin Cleanser



Sempurity DermaPure Antimicrobial Wound and Skin Cleanser is recommended by doctors for the care of chronic wounds. DermaPure wound cleanser is steroid and alcohol free and contains no fragrance, dyes or perfumes. The formula is hypoallergenic and non-irritating to the skin and wound bed. For the management of minor irritations, minor skin abrasions, minor cuts or intact skin.

Item	Product#	Size
Wound Cleanser, Spray Bottle	SNS51118	8oz. btl.

# Wound Care Allowables

## Medicare Allowables...30 day supply

### Medicare Part B

#### Indications and Limitations of Coverage for Surgical Dressings

Alginate Pad .....	30 each
Alginate Pad - with Silver .....	30 each
Alginate Rope (per 6").....	60 inches
Collagen Pad .....	30 each
Collagen Powder .....	30 grams
Foam Pads.....	12 each
Gauze Pad - Sterile with Border.....	30 each
Hydrocolloid Pad .....	12 each
Hydrogel Impregnated Gauze .....	30 each
Hydrogel Pad .....	30 each
Hydrogel Pad with Border.....	12 each
Specialty Absorbent Dressing .....	30 each
Transparent Film Dressing.....	12 each

Safe n Simple has provided the above information as a guideline of Medicare policies and assumes no responsibility for inconsistencies. All medical and claims determinations will be made by Medicare and your insurance plan.

# Wound Care Kits



## Eliminate confusion and ensure your patients receive the products you prescribe

Simpurity wound care kits are your perfect answer - contact us today to see how to:

- Help your patients
- Help your practice
- Help heal wounds

Simpurity wound care kits help you keep control of the dressing choice to achieve optimum wound outcomes, see page 21 for a list of our representatives to call for pricing.

**RECOMMENDED FOR:** Diabetic ulcers, • pressure ulcers, stasis ulcers, post-surgical wounds and donor sites.

**COMPLETE KIT CONTENTS:** Each Kit Includes 15/30 Primary Dressings of your Choice along with 15/30 Simpurity Bordered Gauze Sponge Dressings, 15/30 Packets of Simpurity Gauze Sponges and a bottle of DermaPure Wound and Skin Cleanser.



## KIT CONTAINS:

### PRIMARY DRESSINGS:

**Smpurity Alginate 2"x2" Dressing Pads**  
**Smpurity Silver Alginate 2"x2" Dressing Pads**  
**Smpurity Silver Fibergel 2"x2" Dressing Pads**  
**Smpurity Collagen 2"x2" Dressing Pads**  
**Smpurity Collagen Particles 1g Packets**

Item	Product#	Amount
Alginate 2"x2" Pad Wound Care Kit	SNSWCKA1	15 Day Supply
Alginate 2"x2" Pad Wound Care Kit	SNSWCKA2	30 Day Supply
Silver Alginate 2"x2" Wound Care Kit	SNSWCKS1	15 Day Supply
Silver Alginate 2"x2" Wound Care Kit	SNSWCKS2	30 Day Supply
Silver Fibergel 2"x2" Wound Care Kit	SNSWCKF3	15 Day Supply
Silver Fibergel 2"x2" Wound Care Kit	SNSWCKF4	30 Day Supply
Collagen 2"x2" Wound Care Kit	SNSWCKC1	15 Day Supply
Collagen 2"x2" Wound Care Kit	SNSWCKC2	30 Day Supply
Collagen 1g Pkt Wound Care Kit	SNSWCKP1	15 Day Supply
Collagen 1g Pkt Wound Care Kit	SNSWCKP2	30 Day Supply

# Stay Connected

## Contact Our Team...

---

1. Call Safe n Simple Customer Service  
Toll Free 844-767-6334 or  
Amy Ham 248-214-4877  
Mon-Fri 9am-5pm EST
2. Find a Supplier or Order from our E-Store at:  
[www.sns-medical.com](http://www.sns-medical.com)
3. Billing/Remittance Address:      Shipping Address:  
8112 Rosebud Lane                      5827 Terex  
Clarkston, MI 48348                      Clarkston, MI 48346

## Inside Sales & Ostomy Support Team...

---

Holly Loos                                      248-496-2508  
(AK, AL, AR, AZ, CA, CO, DE, FL, GA, HI, KS, LA, MD, MO, MS, NC, NJ, NM, NV, OK, SC, TN, TX, UT, VA)

Michelle Pitylak                                248-212-2243  
(CT, IA, ID, IL, IN, KY, MA, ME, MI, MN, MT, ND, NE, NH, NY, OH, OR, PA, RI, SD, VA, WA, WI, WV, WY)

## Outside Sales Team...

---

Northeast & Midwest	Kerri Gualtieri	248-802-7084
SE Region Director:	Beth Peter	561-526-6718
Ind. Rep. CA:	Dan Penrose	909-730-0707
Ind. Rep. MD/VA:	Jim Spradlin	804-752-2897
Ind. Rep. Canada:	Peter West	519-726-5959
Bus. Dev./Hospitals:	Tara Franco	248-787-4656

## Samples Team...

---

To request samples call our customer service team  
or email us at: [samples@sns-medical.com](mailto:samples@sns-medical.com)

Clinicians/New Ostomy Patients:  
Request a **SimpleStart Ostomy Sample Kit**  
email us at: [simplestart@sns-medical.com](mailto:simplestart@sns-medical.com) or fax us at: 248-575-4172

# Stay Connected

## Supplier Orders:

---

1. Email Us: [purchaseorder@sns-medical.com](mailto:purchaseorder@sns-medical.com)
2. FAX Us: 248-698-3454
3. Billing/Remittance Address:      Shipping Address:  
8112 Rosebud Lane                      5827 Terex  
Clarkston, MI 48348                      Clarkston, MI 48346

## Marketing Team:

---

Have a special event or sponsorship opportunity you would like to share with Safe n Simple?

Please feel free to email us: [marketing@sns-medical.com](mailto:marketing@sns-medical.com)

Call Us:  
Michelle Bliszack - Marketing Manager  
248-408-9501 Direct  
844-767-6334 Toll Free

**At Safe n Simple we are committed to providing the highest quality products while keeping them affordable to the consumer and providing them the best possible outcomes.**

Log on and “like” our Facebook page.

You will have easy access to promotions, new product updates and much more!

**VISIT US TODAY!**

<https://www.facebook.com/Safe-n-Simple>

Join our Mailing List - Visit: [www.sns-medical.com](http://www.sns-medical.com)

## Our Mission:

Safe n Simple is an innovator, developer and marketer of new medical products. We enhance the patient's quality of life with superior products at an affordable price bringing value to the communities in which we live.

### Safe n Simple Ostomy Care Products

Adhesive Remover  
Peri-Stoma Cleanser & Adhesive Remover  
No-Sting Skin Barrier Film  
Skin Barrier Film  
Ostomy Skin Barrier Powder  
Skin Barrier Ring - Conforming Seal  
Skin Barrier Ring - Extra Thick Conforming Seal  
Skin Barrier Ring - Integrity Seal  
Skin Barrier Ring - Longevity Seal  
Skin Barrier Arc  
X-Tra Wide Skin Barrier Arc  
1/2 Arc Skin Barrier Strips  
Skin Barrier 4"x4" Sheet  
Skin Barrier 4"x4" Sheet w/Aloe & Zinc Oxide  
Skin Barrier Paste  
No-Sting Skin Barrier Paste  
Skin Barrier Paste - Fast Drying  
Ostomy Pouch Deodorant  
Assure C Odor Eliminator  
One-Piece, CTF Pouches  
One-Piece, CTF Pouches with Filters  
Stoma Caps  
Security Hernia/Ostomy Support Belts  
Security Pediatric Bands

### Simpurity Wound Care Products

Alginate - Dressing Pads & Rope  
Alginate Silver - Dressing Pads & Rope  
Collagen - Dressing Pads & Powder  
DermaPro Waterproof Plastic Tape  
DermaPro Waterproof Silicone Tape  
DermaPure Skin and Wound Cleanser  
Fibergel - Dressing Pads  
Fibergel Silver - Dressing Pads  
Foam - Dressing Pads  
Foam - Adhesive Bordered Dressing Pads  
Foam - Fenestrated Bordered Circular Dressing  
Foam - Silicone Adhesive Dressing Pads  
Foam - Silicone Bordered Dressing Pads  
Foam - Silicone Sacrum Bordered Dressing Pads  
Foam - Silicone Silver Foam Dressing Pads  
Foam - Silver Dressing Pads  
Gauzes - Bordered and Non-Bordered  
Hydrocolloid Dressing Pads  
Hydrogel Dressing Sheets  
Hydrogel Gauze Dressing Pads  
IV Derm Transparent Silicone Sheet  
Microroll Gauze  
Polymer Dressing Pads  
Super Absorbent Pads  
Transparent Silicone Film Dressing  
Wound Care Kits  
Xeroform Gauze Dressing Pads



**Safe n Simple**<sup>™</sup>  
Clean. Comfortable. Confident.

Clarkston, MI 48346 • P: 844.767.6334 • F: 248.698.3454

[www.sns-medical.com](http://www.sns-medical.com)

Email: [samples@sns-medical.com](mailto:samples@sns-medical.com)

Like us on Facebook:

<https://www.facebook.com/Safe-n-Simple>

Simpurity<sup>™</sup> is a registered trademark of Safe n Simple, Inc.