



## Product Overview

Safe n Simple focuses on product quality, patient safety, and cost containment. With the introduction of our Collagen wound care powder, we are proud to launch the Simpurity wound care product line. Simpurity wound care products follow evidence based medicine, facilitate moist wound healing and improve patient outcomes while lowering costs.

The Simpurity Collagen in a powder form provides the health care professional with flexibility and adaptability to place collagen in wounds of various sizes and shapes. Simpurity Collagen wound powder is a lyophilized medical grade bovine collagen presenting as a clean, white powder. And, like human connective tissue, Simpurity Collagen wound dressings are rich in Type One Collagen with a unique porous structure that provides excellent absorption and scaffolding for cell proliferation and migration.

Simpurity Collagen provides substrate collagen to divert matrix metalloproteinases (MMP) digestion from healing tissues and reduces inflammation that may slow or prevent healing. Simpurity Collagen fragments provide ingredients for tissue regeneration and the 3D scaffold structure allows fibroblasts to migrate across the wound bed facilitating granulation, angiogenesis and tissue repair. In acute wounds, Simpurity Collagen rapidly activates platelets and coagulating factors to induce hemostasis.



## Product Features

- Easy application and superior conformity to the wound
- Creates a soft, conformable gel on contact with wound fluid
- Highly biocompatible, porous scaffold encourages healing
- Crosslinking process increases mechanical strength and moderates enzymatic digestion
- High exudate absorption and retention provides moist wound healing environment
- May be used as a hemostatic dressing
- 100%, non-bleached, native undigested Bovine Collagen

## Indications

- Partial and full-thickness wounds
- Tunneled / undermined wounds or surgical wounds
- Low to moderately exuding chronic wounds including
  - Diabetic foot ulcers
  - Venous leg ulcers
  - Pressure sores
- Management of light bleeding in acute or surgical wounds
- Post debridement for hemostasis

## Product Details

Product	Item Number	Units/Box	Boxes/Case	HCPC
1 Gram Vial	SNS5001G	5	6	A6010
1 Gram Pack	SNS5221G	10	6	A6010
2"x2" Pad	SNS52222	5	6	A6021

