



Product Overview

Safe n Simple focuses on product quality, patient safety, and cost containment. With the introduction of the Simpurity Foam wound care dressing we are proud to launch the Simpurity wound care product line. The Simpurity wound care products follow evidence based medicine, facilitate moist wound healing and improve patient outcomes while lowering costs.

Simpurity Foam wound dressings are designed to help manage moderate-to-highly draining wounds. The dressing absorbs wound exudate, promotes moist wound healing, protects the wound from contaminants and permits moisture vapor exchange.

The white side of the foam dressing is applied to the wound. The coral colored top of the foam dressing contains a moisture permeable layer designed to protect against exudate strike-through while acting as a protective barrier to outside wound contaminants. Simpurity Foam dressings are designed to be easily removed from the wound sites. The non-adhesive construction of the dressing helps to assure that the skin around the wound perimeter is not damaged.





Product Features

- Easy to remove from wound
- High absorption capacity
- Protective backing to protect against strike through

Indications

Simpurity Foam Wound Dressings are indicated for the management of moderately-to-highly exuding wounds such as:

- Pressure sores
- Leg ulcers
- Cavity wounds
- Use as protective padding

Contra-Indications

Simpurity Foam Wound Dressings are not indicated for use on the following

- Third-degree burns
- Dry wounds
- Surgical implantations

Product Details - Non Adhesive Foam

Product	Item Number	Units/Box	Boxes/Case	HCPC
2" x 2"	SNS51W02	12	16	A6209
4" x 4"	SNS51W04	12	16	A6209
4" x 6"	SNS51W24	12	16	A6210

