

Introducing

# Collagen Pads and Powders

Safe n Simple focuses on product quality, patient safety, and cost containment. Simpurity wound care dressings follow evidence based medicine, facilitates moist wound healing and improves patients outcomes while lowering costs.



**Product Details:**

Size	Item#	Box/Case	HCPC
2"x2"	SNS52222	6 bxs/cs	A6021
4"x4"	SNS52244	6 bxs/cs	A6021
7"x7"	SNS52277	6 bxs/cs	A6023
1g Pkt	SNS5221G	6 bxs/cs	A6010
1g Vial	SNS5001G	6 bxs/cs	A6010

**FEATURES:**

- Triple helix, 100% non-bleached, Type I, native undigested bovine collagen
- Easy application and superior conformity to the wound
- Creates a soft, conformable gel on contact with the wound
- Highly biocompatible, porous scaffold encourages healing
- High exudate absorption and retention provides moist wound healing environment
- May be used as a hemostatic dressing
- Economically priced for a collagen with proven outcomes

**INDICATIONS:**

For the management of low-to-moderately exuding wounds such as:

- Partial and full-thickness wounds
- Tunneled / undermined wounds or surgical wounds
- Low to moderately exuding chronic wounds including:
  - Diabetic foot ulcers
  - Venous leg ulcers
  - Pressure sores
- Management of light bleeding in acute or surgical wounds
- Post debridement for hemostasis

**Create better outcomes and continuity of care for your patients with a quality dressing that has exceptional healing rates!**