



## Product Overview

Safe n Simple focuses on product quality, patient safety, and cost containment. With the introduction of the Simpurity Foam bordered wound care dressing we are proud to expand the Simpurity wound care product line. Simpurity wound care products follow evidence based medicine, facilitate moist wound healing and improve patient outcomes while lowering costs.

Simpurity Foam wound dressings are composed of a soft and a highly absorbent polyurethane foam pad to absorb wound drainage and a transparent adhesive border to hold the dressing in place. This island foam dressing helps maintain a moist environment, which accelerates wound healing. It is especially suitable for moderate to highly exuding wounds.

The white side of the foam dressing with the adhesive is applied to the wound. The clear colored film border of the dressing contains a moisture permeable layer designed to protect against exudate strike-through while acting as a protective barrier to outside wound contaminants. Simpurity Foam dressings adhesive bordered construction helps to assure that the dressing will stay in place and is easily removed from the skin surface surrounding the wound.





## Product Features

- Comfortable and conformable, easy to remove from wound
- High absorption capacity (more than 10 times its weight)
- Transmitting moisture vapor and oxygen and providing moist environment for wound healing
- Semi-permeable polyurethane

## Indications

Simpurity Foam Wound Dressings are indicated for the management of moderately-to-highly exuding wounds such as:

- Pressure sores
- Leg ulcers
- Cavity wounds
- Lacerations and abrasions

## Contra-Indications

Simpurity Foam Wound Dressings are not indicated for use on the following

- Third-degree burns
- Dry wounds
- Surgical implantations

## Product Details - Adhesive Bordered Foam

	;fW °	BSU] SY[ Y	: 5B5
\$ j \$	E@E) \$%" \$	#\$ BSVel 4aj Ž#( 4aj V# 5SeW	3( \$##
& j ' '	E@E) \$%" "	#\$ BSVel 4aj Ž#( 4aj V# 5SeW	3( \$#%
( j (	E@E) \$%" (	#\$ BSVel 4aj Ž#( 4aj V# 5SeW	3( \$#%
* j *	E@E) \$%" *	#\$ BSVel 4aj Ž#( 4aj V# 5SeW	3( \$#&



Clarkston, MI 48348 | 844-767-6334 | [www.SnS-Medical.com](http://www.SnS-Medical.com)

WCF-081022