



Product Overview

Safe n Simple focuses on product quality, patient safety, and cost containment. With the introduction of our Simpurity Silver Alginate wound care dressings we are proud to add to the Simpurity wound care product line. The Simpurity wound care products follow evidence based medicine, facilitate moist wound healing and improve patient outcomes while lowering costs.

Simpurity Silver Alginate dressings are made of natural fibers derived from seaweed through a series of special processes. It is designed to be fast acting and highly absorbant. The dressing absorbs exudates and forms a gel-like covering over the wound, maintaining a moist environment for wound healing. The gel-like surface also keeps the dressing from adhering to the wound so it promotes painless removal. This feature makes dressing changes easier and saves time for both the patient and health care professional. Simpurity Silver Alginate dressing can inhibit bacteria growing in the wound bed and reduce infection. The Simpurity Silver Alginate dressing has been tested in vitro and found to be effective against bacteria such as: Staphylococcus aureus, Enterococcus faecalis, Streptococcus pyogenes, Escherichia coli, Pseudomonas aeruginosa, Klebsiella pneumoniae.



Product Features

- Strong wet strength
- Easy to remove from wound
- High absorption capacity
- Conforms to wound contours
- Reduces Infection
- Gel forming property

Indications

Simpurity Silver Alginate Dressings are indicated for the management of moderately-to-highly exuding partial to full-thickness wounds such as:

- Trauma wounds
- Pressure, Diabetic or Leg ulcers
- Graft and donor sites
- Post operative surgical wounds

Contra-Indications

Simpurity Silver Alginate Dressings are not indicated for use on the following:

- Patients with known sensitivity to silver
- Third degree burns

Technical Specifications

| | |
|------------------------------|--|
| Antifungal Agent: | Silver Particles |
| Active against pathogens: | Gram positive and gram negative bacteria |
| Spectrum of activity: | Broad spectrum antibacterial effect for gram positive and gram negative bacteria |
| Concentration on the device: | Each dressing contains 50ug/100cm ² silver particles. |

Product Details

| Size | Item# | Packaging | HCPC |
|-----------------|----------|----------------------------------|-------|
| 2" x 2" | SNS51702 | 10 Dressings/Box - 16 Boxes/Case | A6196 |
| 4" x 4" | SNS51704 | 10 Dressings/Box - 16 Boxes/Case | A6196 |
| 4" x 5" | SNS51720 | 10 Dressings/Box - 16 Boxes/Case | A6197 |
| 4" x 6" | SNS51724 | 10 Dressings/Box - 16 Boxes/Case | A6197 |
| 6" x 6" | SNS51706 | 10 Dressings/Box - 16 Boxes/Case | A6197 |
| 8" x 8" | SNS51764 | 10 Dressings/Box - 16 Boxes/Case | A6198 |
| 1" x 12" (rope) | SNS51712 | 10 Dressings/Box - 12 Boxes/Case | A6199 |



Clarkston, MI 48348 | 844-767-6334 | www.SnS-Medical.com