
Introducing

Fibergel

with Gelling Action

Safe n Simple focuses on product quality, patient safety, and cost containment. Simpurity wound care dressings follow evidence based medicine, facilitates moist wound healing and improves patients outcomes while lowering costs.



Product Details:

Size	Item#	Box/Case	HCPC
2"x2"	SNS56702	10 bx/16 bxs	A6196
4"x4"	SNS56704	10 bx/16 bxs	A6196
4"x5"	SNS56745	10 bx/16 bxs	A6197

FEATURES:

- Soft, absorbent material that forms into a gel when in contact with exudate
- Gelling action creates an optimal environment for moist wound healing
- Conforms easily to various wound shapes
- Supports debridement
- Can easily be removed in one piece with minimal damage to newly forming granulating tissue
- For use as a primary sterile dressing

INDICATIONS:

For the management of moderate-to-highly exuding wounds such as:

- Pressure Ulcers
- Leg Ulcers
- Cavity Wounds
- Post Debridement Hemostasis
- Superficial and Partial Thickness Wounds
- Traumatic Wounds

Comparison testing has demonstrated that Simpurity Fibergel wound dressing has an excellent performance against leading hydrofiber dressings in:

Absorbency - absorbs a high level of exudate during usage **Fluid Retention Capacity -** retains absorbed fluids under compression **Wet Strength -** remains intact when being removed from the wound bed

