

*Introducing*

# HydroGel

## Absorbent Sheet Wound Dressings

Safe n Simple focuses on product quality, patient safety, and cost containment. Simpurity wound care dressings follow evidence based medicine, facilitates moist wound healing and improves patients outcomes while lowering costs.



**FEATURES:**

- Provides a moist wound healing environment for wounds with minimal to no exudate
- Permeable to air and moisture vapor for faster wound healing
- Sterile and Latex Free
- Gentle to remove without damaging tissues
- Provides a soothing coolness to the wound bed
- Can be cut to fit the wound site
- Conforms easily to the wound bed
- Prevents bacterial infection and contamination

**INDICATIONS:**

For the management of wounds such as:

- Skin Tears
- Superficial abrasions
- 1st and 2nd degree burns
- Stage I or II pressure ulcers
- Radiation tissue damage

**Product Details:**

**HydroGel Dressing Sheets - Non Bordered**

Size	Item#	Box	HCPC
2"x2"	SNS58302	10 Sheets	A6242
4"x5"	SNS58345	10 Sheets	A6243

**HydroGel Dressing Sheets - Bordered**

Size	Item#	Box	HCPC
2"x2"	SNS58312	12 Sheets	A6245
4"x5"	SNS58315	12 Sheets	A6246

**RECOMMENDED FOR:**

Acute to chronic wound care