



Wound Care Product Guide

Introduction...

Safe n Simple™ is a medical care innovator, specializing in the promotion of healthy skin, complete security and confident ostomy and wound care management. We view our product development through the eyes of our customers, who want active, productive lives unimpeded by their conditions. Our goal is to make our customer's ostomy and wound care safe and simple.

We follow two overriding principles: create innovative, high-performance products, and deliver the greatest possible value. We continually push our technicians to develop better ideas, more advanced formulations and superior products, while at the same time making those products convenient and affordable.

We've spent years in **Ostomy care** learning from our customers what they really want. The result is a complete line of ostomy care products designed to *Cleanse, Protect, Secure, and Assure*.

And now a complete line of wound care products created to promote healing and cost savings to users.

Above all, our customers want safety and simplicity. With Safe n Simple, that's our promise.

Table of Contents

Wound Care

Page 2-4	Wound Care Intro Pages & Table of Contents
Page 5	Alginate
.....	Silver Alginate
Page 6	Collagen
.....	Fibergel
Page 7	Silver Fibergel
.....	Borderless Foam
Page 8	Bordered Foam
.....	Ionic Silver Foam
Page 9	Silicone Sacrum Bordered Foam
.....	Borderless Silicone Foam
Page 10.....	Bordered Silicone Foam
.....	Borderless Silver Silicone Foam
Page 11.....	Bordered Silver Silicone Foam
.....	Gauze Sponges
Page 12.....	Bordered Gauze Sponges
.....	Hydrocolloids
Page 13.....	Hydrogel Bordered Sheets
.....	Hydrogel Sheets
Page 14.....	Hydrogel Impregnated Gauze
.....	Polymer Pads
Page 15.....	Silicone Transparent Film
.....	IV Derm Silicone Transparent Film
Page 16.....	DermaPro Waterproof Plastic Tape
.....	DermaPro Waterproof Silicone Tape
Page 17.....	Xeroform Impregnated Gauze
.....	DermaPure Antimicrobial Wound and Skin Cleanser
Page 18 & 19.....	Wound Care Kits
Page 20.....	Wound Care Allowables
Page 21-23	Stay Connected with SNS



Safe n Simple WOUND CARE
Clean. Comfortable. Confident.



Enhancing quality of life with exceptional products at an affordable price

Safe n Simple focuses on product quality, patient safety, and cost containment. With the introduction of the Simpurity Wound Care product line we are proud to offer wound care products that follow evidence based medicine, facilitate moist wound healing and improve patient outcomes while lowering costs.

Your patients deserve the best care available, with the proven quality of Simpurity wound care product solutions, you can help them achieve greater outcomes and save money doing it.

Testimonial.....

"I personally use Safe n Simple wound care products in our wound center for several reasons. The Simpurity line of products offers high quality alginate, foam and collagen dressings with and without silver at a very cost effective price. When everyone is looking to use evidenced based treatments and work within our current reimbursement environment this is the PERFECT product line. The products are user friendly. The packaging is appealing and instructions clear. Most importantly our healing rates and results have been excellent."

J. Christopher Brandys, MD FRCS FACS
*General and Vascular Surgeon
Director of The Wound Center
Medical Consultant - Safe n Simple*

For more information or to request samples please call our toll-free number or visit our website.

844-767-6334 | www.SnS-Medical.com

Affordable products covered by insurance.

Alginates

Alginate Wound Dressing

Simpurity Alginate Dressings are derived from seaweed. The products are composed of soft, calcium alginate fibers shaped as pads (fibrous mats) or ropes (twisted fibers). When packed in a wound, the alginate interacts with wound exudate to form a soft gel that maintains a moist healing environment. The dressings release of calcium ions can activate prothrombin and improve hemostasis. Has a strong tensile strength.



Item	Product#	Size	HCPC
Dressing Pad, 10/Box	SNS50702	2"x2"	A6196
Dressing Pad, 10/Box	SNS50704	4"x4"	A6196
Dressing Pad, 10/Box	SNS50720	4"x5"	A6197
Dressing Pad, 5/Box	SNS50732	4"x8"	A6197
Dressing Pad, 10/Box	SNS50764	8"x8"	A6198
Rope, 30/Bag	SNS50712	12"	A6199

Silver Alginate Wound Dressing

Simpurity Silver Alginate Dressings are made of natural fibers derived from seaweed through a series of special processes. Simpurity Silver Alginate Dressings can inhibit bacteria growing in the wound bed reducing risk of infection while maintaining a moist environment for wound healing. The dressing has a strong tensile strength, a high absorption capacity and is easy to remove from the wound bed.



Item	Product#	Size	HCPC
Dressing Pad, 10/Box	SNS51702	2"x2"	A6196
Dressing Pad, 10/Box	SNS51704	4"x4"	A6196
Dressing Pad, 10/Box	SNS51720	4"x5"	A6197
Dressing Pad, 10/Box	SNS51764	8"x8"	A6198
Rope, 10/Box	SNS51712	12"	A6199

Collagens

Collagen Wound Dressing

Simpurity Collagen provide ingredients for tissue regeneration in phase 3 wounds and the 3D scaffold structure of Simpurity Collagen allows fibroblasts to migrate across the wound bed facilitating granulation, angiogenesis and tissue repair. Simpurity Collagen dressings are 100% pure Type I, non-bleached, native undigested bovine collagen. The dressings provide high exudate absorption and retention providing a moist wound healing environment. Easy application and superior conformity to the wound.



Item	Product#	Size	HCPC
Dressing Pad, 10/Box	SNS50002	2"x2"	A6021
Dressing Pad, 5/Box	SNS52222	2"x2"	A6021
Dressing Pad, 5/Box	SNS52244	4"x4"	A6021
Vial Particles, 5/Box	SNS5001G	1 Gram	A6010
Packet Particles, 10/Box	SNS5221G	1 Gram	A6010

Note: SNS50002 is a premium brand collagen wound dressing.

Fibergels

Fibergel Wound Dressing

Simpurity Fibergel Wound Dressing is a soft, absorbent material that forms into a gel when in contact with exudate. The gelling action creates an optimal environment for moist wound healing. It conforms easily to various wound shapes, supports debridement and can easily be removed in one piece with minimal damage to newly forming granulating tissues.



Item	Product#	Size	HCPC
Dressing Pad, 10/Box	SNS56702	2"x2"	A6196
Dressing Pad, 10/Box	SNS56704	4"x4"	A6196
Dressing Pad, 10/Box	SNS56745	4"x5"	A6197

Fibergels

Silver Fibergel Wound Dressing

Simpurity Fibergel Ag Wound Dressing has all the benefits of a hydrofiber dressing plus ionic silver to inhibit bacteria from growing in the wound bed, reducing the risk of infection. The gelling action of Fibergel AG creates an optimal environment for moist wound healing. It conforms easily to various wound shapes, supports debridement and can easily be removed in one piece with minimal damage to newly forming granulating tissues.



Item	Product#	Size	HCPC
Dressing Pad, 10/Box	SNS56712	2"x2"	A6196
Dressing Pad, 10/Box	SNS56744	4"x4"	A6196
Dressing Pad, 10/Box	SNS56716	4"x4.75"	A6197

Foams

Foam Non-Adhesive Wound Dressing

Simpurity Foam Wound Dressings are designed to help manage moderate-to-highly draining wounds. The dressing absorbs wound exudate, promotes moist wound healing, protects the wound from contaminants and permits moisture vapor exchange. These foam dressings have a protective colored film backing to protect against strike through.



Item	Product#	Size	HCPC
Dressing Pad, 12/Box	SNS51W02	2"x2"	A6209
Dressing Pad, 12/Box	SNS51W03	3"x3"	A6209
Dressing Pad, 12/Box	SNS51W04	4"x4"	A6209
Dressing Pad, 12/Box	SNS51W24	4"x6"	A6210

Foams

Foam Wound Dressing - Adhesive Border

Simpurity Foam wound dressings are composed of a soft and highly absorbent polyurethane foam pad to absorb wound drainage and a transparent adhesive border to hold the dressing in place. This island foam dressing helps maintain a moist environment, which accelerates wound healing. It is especially suitable for moderate to highly exuding wounds and has a protective backing to protect against strike through. The dressings are comfortable, easy to apply, and easy to remove from the wound.



Item	Product#	Size	HCPC
Dressing Pad, 12/Box	SNS72302	2"x2"	A6212
Dressing Pad, 12/Box	SNS72320	4"x5"	A6213

Foam Wound Dressing - Ionic Silver

Simpurity Foam Dressing is a silver ion-containing polyurethane foam dressing that protects wounds from bacterial penetration while reducing the bioburden in a wound bed. It mainly consists of hydrophilic polyurethane foam substrate and silver compound. It is indicated for the management of moderate-to-high exudate for chronic and acute wounds and helps maintain a moist healing environment. Controlled release of ionic silver provides sustained antibacterial effect. Inhibits a broad-spectrum of infection up to 99.99% for a wide variety of bacteria. Provides longer wear times and pain-free removal.



Item	Product#	Size	HCPC
Dressing Pad, 12/Box	SNS72502	2"x2"	A6209
Dressing Pad, 12/Box	SNS72520	4"x5"	A6210

Foams

Silicone Sacrum Foam Wound Dressing

Simpurity Sacrum Foam wound dressings are composed of a soft absorbent pad that enhances patient comfort and absorbs excess fluids. The gentle silicone border adheres to the surrounding skin, not the wound bed. Helps to minimize pain and trauma at dressing changes. Minimizes pressure ulcers of the sacral region.



Item	Product#	Size	HCPC
Dressing Pad, 12/Box	SNS77777	7"x7"	A6213

Foam Wound Dressing - Silicone

Simpurity Silicone Foam Dressing is a superior absorbent dressing for moderate to high wound exudate management. It contains a highly absorbent foam pad and gentle adhesive silicone contact layer. These dressings help maintains a moist wound healing environment with no irritation upon removal, minimizes trauma to the wound and reduces pain during dressing changes.



Item	Product#	Size	HCPC
Dressing Pad, 12/Box	SNS74502	2"x2"	A6209
Dressing Pad, 12/Box	SNS74520	4"x5"	A6210

Foams

Foam Wound Dressing - Bordered Silicone

Simpurity Silicone Foam Bordered Dressing is a superior absorbent dressing for moderate to high wound exudate management. It contains a highly absorbent foam pad and gentle adhesive silicone contact layer and an adhesive border to hold the dressing in place. These dressings help maintain a moist wound healing environment with no irritation upon removal, minimizes trauma to the wound and reduces pain during dressing changes.



Item	Product#	Size	HCPC
Dressing Pad, 12/Box	SNS77733	3"x3"	A6212
Dressing Pad, 12/Box	SNS77744	4"x4"	A6212

Foam Wound Dressing - Silver Silicone

Simpurity Silicone Foam Dressing is a polyurethane foam pad containing silver ions that has a gentle adhesive silicone contact layer to protect wounds from bacterial penetration while reducing the bioburden in the wound bed. It is indicated for the management of moderate-to-high exudate for chronic and acute wounds and helps maintain a moist healing environment. Controlled release of ionic silver provides sustained antibacterial effect. Provides longer wear times and pain-free removal.



Item	Product#	Size	HCPC
Dressing Pad, 12/Box	SNS74422	2"x2"	A6209
Dressing Pad, 12/Box	SNS74445	4"x5"	A6210
Dressing Pad, 12/Box	SNS74466	6"x6"	A6210
Dressing Pad, 12/Box	SNS74477	7"x7"	A6211

Foams

Foam Wound Dressing - Bordered Silver Silicone

Simpurity Silver Silicone Foam Bordered Dressing is a superior absorbent dressing containing silver ions for moderate to high wound exudate management. It contains a highly absorbent foam pad impregnated with silver ions and gentle adhesive silicone contact layer with an adhesive border to hold the dressing in place. These dressings help maintain a moist wound healing environment with no irritation upon removal, minimizes trauma to the wound and reduces pain during dressing changes. The 2"x2" size has a 1" sq. pad, the 4.6"x4.6" has a 3" sq. pad.



Item	Product#	Size	HCPC
Dressing Pad, 12/Box	SNS77111	2"x2"	A6212
Dressing Pad, 12/Box	SNS77133	4.6"x4.6"	A6212

Gauzes

Gauze Sterile Sponge

Simpurity Gauze Sterile Sponges are ideal for general wound cleansing, debriding, minor prepping, covering and packing wounds. It is a woven sponge made of 12 ply 100% USP cotton. Contains 2 sponges per packet.



Item	Product#	Size	HCPC
Dressing Pad, 2 Pcs/pkt, 25/Box	SNS54444	4"x4"	A6402

Gauzes

Bordered Gauze Sterile Sponge

Simpurity Bordered Gauze Sterile Sponges are a non-adherent gauze sponge island dressing with an adhesive border. These woven sponges are made of 12ply 100% USP cotton for use as a primary or secondary dressing. The sponge absorbs exudate for a cleaner wound bed and the adhesive border protects the wound from outside contaminants. The 4"x4" has gauze pad size of 2"x2", and 6"x6.75" has gauze pad size of 4"x4.75".



Item	Product#	Size	HCPC
Dressing Pad, 30/Box	SNS52316	4"x4"	A6219
Dressing Pads, 30/Box	SNS52324	6"x6.75"	A6220

Hydrocolloids

Hydrocolloid Wound Dressing

Simpurity Hydrocolloid Wound Dressing is a thin hydrocolloid sheet. Simpurity Hydrocolloids are absorbent, and delivers a gentle and secure adhesion that provides a barrier to external debris and bacteria when securely adhered to the skin. For use on light-to-moderately exuding wounds.



Item	Product#	Size	HCPC
Dressing Pad, 12/Box	SNS55402	2"x2"	A6234
Dressing Pad, 12/Box	SNS55404	4"x4"	A6234
Dressing Pad, 10/Box	SNS21149	7"x7"	A6236

Hydrogels

Hydrogel Bordered Dressing Sheet


Simpurity Hydrogel Absorbent Dressing Sheets have a polyurethane border and are designed to provide a moist wound healing environment for wounds with minimal to no exudate. It is permeable to air and moisture vapor for faster wound healing. Easy to remove from the skin without damaging tissues. Conforms to the wound bed preventing bacterial infection and contamination.



Item	Product#	Size	HCPC
Dressing Pad, 12/Box	SNS58312	2"x2"	A6245
Dressing Pad, 12/Box	SNS58315	4"x5"	A6246

Hydrogel Dressing Sheet

Simpurity Hydrogel Absorbent Dressing Sheets are designed to provide a moist wound healing environment for wounds with minimal to no exudate. It is permeable to air and moisture vapor for faster wound healing. Easy to remove from the skin without damaging tissues. Conforms to the wound bed preventing bacterial infection and contamination.



Item	Product#	Size	HCPC
Dressing Pad, 10/Box	SNS58302	2"x2"	A6242
Dressing Pad, 10/Box	SNS58345	4"x5"	A6243

Hydrogels

Hydrogel - Impregnated Gauze

Simpurity Hydrogel Impregnated Gauze Wound Dressing is designed to provide a moist wound healing environment for wounds with minimal or no exudate. Contains a high viscosity gel impregnated into the gauze pad. Easy to apply.



Item	Product#	Size	HCPC
Dressing Pad, 10/Box	SNS58802	2"x2"	A6231
Dressing Pad, 10/Box	SNS58820	4"x5"	A6232

Polymer Pads

Polymer Wound Dressing - Highly Absorbent

Simpurity Polymer Foam Wound Dressing is a thin, sterile, superabsorbent single-use dressing, for use on highly exuding wounds. Simpurity Polymer Foam absorbs the exudate and locks it into the super-absorbent core, allowing for significantly more fluid to be handled by the dressing while keeping it away from the wound and surrounding tissue. This can help to reduce maceration, as well as reduce exposure to harmful components of wound fluid.



Item	Product#	Size	HCPC
Dressing Pad, 10/Box	SNS59020	4"x5"	A6251
Dressing Pad, 10/Box	SNS59080	8"x10"	A6253

Silicone Transparent Films

Transparent Silicone Film Wound Dressing

Simpurity Transparent Silicone Film wound dressings provide a waterproof barrier to external contaminants. It is permeable to air and moisture vapor for faster wound healing. Gentle to remove from sensitive skin areas, can be changed without pain or damage to the skin.



Item	Product#	Size	HCPC
Dressing Sheet, 12/Box	SNS57245	4"x5"	A6258

IV Derm Silicone Film Dressing

Simpurity IV Derm Transparent Silicone Film dressings provide a waterproof barrier to secure IV lines from outside contaminants. Gentle to remove from sensitive skin areas, can be changed without pain or damage to the skin.



Item	Product#	Size	HCPC
Dressing Sheet, 12/Box	SNS52723	2"x3"	A6257

Tapes

DermaPro Waterproof Plastic Tape

Simpurity DermaPro Waterproof Plastic Tape is superior for securing wound dressings and repelling moisture from outside contaminants. DermaPro has a latex free adhesive and is gentle to remove from sensitive skin areas. The roll contains a start strip for easy dispensing.



Item	Product#	Size	HCPC
Plastic Tape, 1 Roll/Box	SNS57220	2"x15'	A4452

DermaPro Waterproof Silicone Tape

Simpurity DermaPro Waterproof Silicone Tape is latex free and has a silicone adhesive that is gentle to remove from sensitive skin areas. Superior for securing wound dressings and repelling moisture from outside contaminants. Each roll is perforated for easy tearing. Contains 12-1" rolls per box or 6-2" rolls per box.



Item	Product#	Size	HCPC
Silicone Tape, 12 Rolls/Box	SNS57230	1"x15'	A4452
Silicone Tape, 6 Rolls/Box	SNS57232	2"x15'	A4452

Xeroforms

Xeroform - Impregnated Gauze

Simpurity Xeroform Impregnated Gauze is a sterile, occlusive, non-adhering dressing consisting of an absorbent, fine-mesh gauze impregnated with an antimicrobial petrolatum blend. The dressing will conform to all body contours which makes it suitable for use in compression therapy, deep and shallow wounds, burns, plastic surgery and surgical wounds with light exudate where mild medication and deodorization are desired.



Item	Product#	Size	HCPC
Dressing Pad, 10/Box	SNS58818	1"x8"	A6222
Dressing Pad, 10/Box	SNS58822	2"x2"	A6222
Dressing Pad, 10/Box	SNS58844	4"x4"	A6222
Dressing Pad, 10/Box	SNS58859	5"x9"	A6223

Wound and Skin Cleanser

DermaPure Antimicrobial Wound and Skin Cleanser



Simpurity DermaPure Antimicrobial Wound and Skin Cleanser is recommended by doctors for the care of chronic wounds. DermaPure wound cleanser is steroid and alcohol free and contains no fragrance, dyes or perfumes. The formula is hypoallergenic and non-irritating to the skin and wound bed. For the management of minor irritations, minor skin abrasions, minor cuts or intact skin.

Item	Product#	Size
Wound Cleanser, Spray Bottle	SNS51118	8oz. btl.

Wound Care Kits



Eliminate confusion and ensure your patients receive the products you prescribe

Simpurity wound care kits are your perfect answer - contact us today to see how to:

- Help your patients
- Help your practice
- Help heal wounds

Simpurity wound care kits help you keep control of the dressing choice to achieve optimum wound outcomes, see page 21 for a list of our representatives to call for pricing.

RECOMMENDED FOR: Diabetic ulcers, • pressure ulcers, stasis ulcers, post-surgical wounds and donor sites.

COMPLETE KIT CONTENTS: Each Kit Includes 15/30 Primary Dressings of your Choice along with 15/30 Simpurity Bordered Gauze Sponge Dressings, 15/30 Packets of Simpurity Gauze Sponges and a bottle of DermaPure Wound and Skin Cleanser.



KIT CONTAINS:

PRIMARY DRESSINGS:

- Simpurity Alginate 2"x2" Dressing Pads
- Simpurity Silver Alginate 2"x2" Dressing Pads
- Simpurity Silver Fibergel 2"x2" Dressing Pads
- Simpurity Collagen 2"x2" Dressing Pads
- Simpurity Collagen Particles 1g Packets

Item	Product#	Amount
Alginate 2"x2" Pad Wound Care Kit	SNSWCKA1	15 Day Supply
Alginate 2"x2" Pad Wound Care Kit	SNSWCKA2	30 Day Supply
Silver Alginate 2"x2" Wound Care Kit	SNSWCKS1	15 Day Supply
Silver Alginate 2"x2" Wound Care Kit	SNSWCKS2	30 Day Supply
Silver Fibergel 2"x2" Wound Care Kit	SNSWCKF3	15 Day Supply
Silver Fibergel 2"x2" Wound Care Kit	SNSWCKF4	30 Day Supply
Collagen 2"x2" Wound Care Kit	SNSWCKC1	15 Day Supply
Collagen 2"x2" Wound Care Kit	SNSWCKC2	30 Day Supply
Collagen 1g Pkt Wound Care Kit	SNSWCKP1	15 Day Supply
Collagen 1g Pkt Wound Care Kit	SNSWCKP2	30 Day Supply

Wound Care Allowables

Medicare Allowables...30 day supply

Medicare Part B - Indications and Limitations of Coverage for Surgical Dressings

Alginate Pad	30 each
Alginate Pad - with Silver.....	30 each
Alginate Rope (per 6").....	60 inches
Collagen Pad	30 each
Collagen Powder	30 grams
Foam Pads.....	12 each
Gauze Pad - Sterile with Border.....	30 each
Hydrocolloid Pad	12 each
Hydrogel Impregnated Gauze	30 each
Hydrogel Pad	30 each
Hydrogel Pad with Border	12 each
Specialty Absorbent Dressing	30 each
Transparent Film Dressing.....	12 each

Safe n Simple has provided the above information as a guideline of Medicare policies and assumes no responsibility for inconsistencies. All medical and claims determinations will be made by Medicare and your insurance plan.

Stay Connected

Contact Our Team...

1. Call Safe n Simple Customer Service
Toll Free 844-767-6334 or
Amy Ham 248-214-4877
Mon-Fri 9am-5pm EST
2. Find a Supplier or Order from our E-Store at:
www.sns-medical.com
3. Billing/Remittance Address: Shipping Address:
8112 Rosebud Lane 5827 Terex
Clarkston, MI 48348 Clarkston, MI 48346

Inside Sales & Ostomy Support Team...

Holly Loos 248-496-2508
(AK, AL, AR, AZ, CA, CO, DE, FL, GA, HI, KS, LA, MD, MO, MS, NC, NJ, NM, NV, OK, SC, TN, TX, UT, VA)

Michelle Pitylak 248-212-2243
(CT, IA, ID, IL, IN, KY, MA, ME, MI, MN, MT, ND, NE, NH, NY, OH, OR, PA, RI, SD, VA, WA, WI, WV, WY)

Outside Sales Team...

Northeast & Mideast	Kerri Gualtieri	248-802-7084
SE Region Director:	Beth Peter	561-526-6718
AR, KY, TN, South TX:	Joel Holt	615-739-3435
Ind. Rep. CA:	Dan Penrose	909-730-0707
Ind. Rep. MD/VA:	Jim Spradlin	804-752-2897
Ind. Rep. Canada:	Peter West	519-726-5959
GPO's/Hospitals:	Tara Franco	248-787-4656

Samples Team...

To request samples call our customer service team
or email us at: samples@sns-medical.com

Clinicians/New Ostomy Patients:
Request a **SimpleStart Ostomy Sample Kit**
email us at: simplestart@sns-medical.com or fax us at: 248-575-4172

Stay Connected

Supplier Orders:

1. Email Us: purchaseorder@sns-medical.com
2. FAX Us: 248-698-3454
3. Billing/Remittance Address: Shipping Address:
8112 Rosebud Lane 5827 Terex
Clarkston, MI 48348 Clarkston, MI 48346

Marketing Team:

Have a special event or sponsorship opportunity you would like to share with Safe n Simple?

Please feel free to email us: marketing@sns-medical.com

Call Us:
Michelle Bliszack - Marketing Manager
248-408-9501 Direct
844-767-6334 Toll Free

At Safe n Simple we are committed to providing the highest quality products while keeping them affordable to the consumer and providing them the best possible outcomes.

Log on and “like” our Facebook page.

You will have easy access to promotions, new product updates and much more!

VISIT US TODAY!

<https://www.facebook.com/Safe-n-Simple>

Join our Mailing List - Visit: www.sns-medical.com

Safe n Simple
Clean. Comfortable. Confident.



Our Mission:

Safe n Simple is an innovator, developer and marketer of new medical products. We enhance the patient's quality of life with superior products at an affordable price bringing value to the communities in which we live.

Safe n Simple Ostomy Care Products

Adhesive Remover
Peri-Stoma Cleanser & Adhesive Remover
No-Sting Skin Barrier Film
Skin Barrier Film
Ostomy Skin Barrier Powder
Skin Barrier Ring - Conforming Seal
Skin Barrier Ring - Integrity Seal
Skin Barrier Ring - Longevity Seal
Skin Barrier Arc
X-Tra Wide Skin Barrier Arc
1/2 Arc Skin Barrier Strips
Skin Barrier 4"x4" Sheet
Skin Barrier Paste
No-Sting Skin Barrier Paste
Skin Barrier Paste - Fast Drying
Ostomy Pouch Deodorant
Assure C Odor Eliminator
One-Piece, CTF Pouches
One-Piece, CTF Pouches with Filters
Stoma Caps
Security Hernia/Ostomy Support Belts
Security Pediatric Bands

Simpurity Wound Care Products

Alginate - Dressing Pads & Rope
Alginate Silver - Dressing Pads & Rope
Collagen - Dressing Pads & Powder
DermaPro Waterproof Plastic Tape
DermaPro Waterproof Silicone Tape
DermaPure Skin and Wound Cleanser
Fibergel - Dressing Pads
Fibergel Silver - Dressing Pads
Foam - Dressing Pads
Foam - Adhesive Bordered Dressing Pads
Foam - Silicone Adhesive Dressing Pads
Foam - Silicone Bordered Dressing Pads
Foam - Silicone Sacrum Bordered Dressing Pads
Foam - Silicone Silver Foam Dressing Pads
Foam - Silver Dressing Pads
Hydrocolloid Dressing Pads
Hydrogel Dressing Sheets
Hydrogel Gauze Dressing Pads
IV Derm Transparent Silicone Sheet
Polymer Dressing Pads
Transparent Silicone Film Dressing
Wound Care Kits
Xeroform Gauze Dressing Pads



Safe n Simple[™]
Clean. Comfortable. Confident.

Clarkston, MI 48346 • P: 844.767.6334 • F: 248.698.3454

www.sns-medical.com

Email: samples@sns-medical.com

Like us on Facebook:

<https://www.facebook.com/Safe-n-Simple>

Simpurity[™] is a registered trademark of Safe n Simple, Inc.