



# **Product Overview**

Safe n Simple focuses on product quality, patient safety, and cost containment. With the introduction of the Simpurity Polymer Pad wound care dressing we are proud to expand the Simpurity wound care product line. Simpurity wound care products follow evidence based medicine, facilitate moist wound healing and improve patient outcomes while lowering costs.

Simpurity Polymer Pad wound dressings are composed of a soft and highly absorbent polymer pad, with a white mesh contact layer and light blue adhesive free outer cloth protection layer. This polymer pad wound dressing helps maintain a moist environment by absorbing proteases and MMPs while wicking excess wound drainage and harmful components away from the wound and traps them in the dressing to accelerate healing. It is especially suitable for moderate to highly exuding wounds.

The white mesh side of the polymer pad wound dressing is applied to the wound. The dressing contains a moisture permeable adhesive free blue outer cloth layer designed to protect against exudate strike-through while acting as a protective barrier to outside wound contaminants. Simpurity Polymer Pad wound dressings bordered construction helps to assure that the dressing will be comfortable, stay in place, and is easily removed from the wound surface.





#### **Product Features**

- Comfortable and conformable, easy to remove from wound
- High absorption capacity (more than 10x its weight)
- · Helps control maceration & excoriation of the wound edge to promote healing
- · Absorbs exudate to help reduce odor and toxins from the wound
- Can be used under compression

# Indications

Simpurity Polymer Pad Wound Dressings are indicated for the management of moderate-to-highly exuding wounds such as:

- Pressure sores
- Venous ulcers
- First & second degree burns
- Post-operative wounds

### **Contra-Indications**

Simpurity Polymer Pad Wound Dressings are not indicated for use on the following

- Third-degree burns
- Dry wounds
- Surgical implantations

#### **Product Details** - Polymer Pad Wound Dressing

Product	Item Number	Units/Box	Boxes/Case	НСРС
4" × 5"	SNS59020	10	8	A6251
8" × 10"	SNS59080	10	8	A6253



Clarkston, MI 48346 | 844-767-6334 | www.SnS-Medical.com