



## Product Overview

Safe n Simple focuses on product quality, patient safety, and cost containment. With the introduction of our Simpurity Silver Alginate wound care dressings we are proud to add to the Simpurity wound care product line. The Simpurity wound care products follow evidence based medicine, facilitate moist wound healing and improve patient outcomes while lowering costs.

Simpurity Silver Alginate dressings are made of natural fibers derived from seaweed through a series of special processes. It is designed to be fast acting and highly absorbant. The dressing absorbs exudates and forms a gel-like covering over the wound, maintaining a moist environment for wound healing. The gel-like surface also keeps the dressing from adhering to the wound so it promotes painless removal. This feature makes dressing changes easier and saves time for both the patient and health care professional. Simpurity Silver Alginate dressing can inhibit bacteria growing in the wound bed and reduce infection. The Simpurity Silver Alginate dressing has been tested in vitro and found to be effective against bacteria such as: Staphylococcus aureus, Enterococcus faecalis, Streptococcus pyogenes, Escherichia coli, Pseudomonas aeruginosa, Klebsiella pneumoniae. Simpurity Silver Alginate dressings are available in 2"x2", 4"x4", 4"x5" and 8"x8" pad sizes, and 1"x12" Rope.



## Product Features

- Strong wet strength
- Easy to remove from wound
- High absorption capacity
- Conforms to wound contours
- Reduces Infection
- Gel forming property

## Indications

Simpurity Silver Alginate Dressings are indicated for the management of moderately-to-highly exuding partial to full-thickness wounds such as:

- Trauma wounds
- Pressure, Diabetic or Leg ulcers
- Graft and donor sites
- Post operative surgical wounds

## Contra-Indications

Simpurity Silver Alginate Dressings are not indicated for use on the following:

- Patients with known sensitivity to silver
- Third degree burns

## Technical Specifications

Antifungal Agent:	Silver Particles
Active against pathogens:	Gram positive and gram negative bacteria
Spectrum of activity:	Broad spectrum antibacterial effect for gram positive and gram negative bacteria
Concentration on the device:	Each dressing contains 50ug/100cm <sup>2</sup> silver particles.

## Product Details

Product	Item Number	Units/Box	Boxes/Case	HCPC
2" x 2"	SNS51702	10	16	A6196
4" x 4"	SNS51704	10	16	A6196
4" x 5"	SNS51720	10	16	A6197
8" x 8"	SNS51764	10	16	A6198
1" x 12" Rope	SNS51712	10	12	A6199



Clarkston, MI 48348 | 844-767-6334 | [www.SnS-Medical.com](http://www.SnS-Medical.com)

Simpurity™ is a registered trademark of Safe n Simple, LLC.