

Introducing

Microroll™ Gauze Wound Dressing

Safe n Simple focuses on product quality, patient safety, and cost containment. Simpurity wound care dressings follow evidence based medicine, facilitates moist wound healing and improves patients outcomes while lowering costs.



Product Details:

Size: 4.5"x4.1yds	Item: SNS53344	Case: 100 Rolls	HPCP: A6266
-----------------------------	--------------------------	---------------------------	-----------------------

FEATURES:

- Absorbs exudate - fast wicking action
- Superior breathability
- Low cost
- Contains PHMB as an antimicrobial agent for inhibiting bacterial growth
- Attacks bacteria in and on the dressing
- Provides prophylactic infection control
- Hypoallergenic, sterile and latex free

INDICATIONS:

For use as a primary or secondary dressing for the management of wounds such as:

- Acute and chronic wounds
- Moderate to highly exuding wounds
- Deep wounds
- Cavity wounds
- Burns

In this study the wounds of five subjects were treated with five daily changes of the PHMB dressing and a control group of five was treated with a non-PHMB dressing. Wounds were swabbed daily and the presence or absence of four specific wound pathogens were determined. The results indicated a marked decrease in frequency of isolation of MRSA and P. aeruginosa when treated with the drain sponge containing PHMB.

Maximum number of days possible for each organism = 25 (5 patients tested 5 times). From Motta and Triglia, 2005

Pathogen	Number of days pathogen present	
	PHMB dressing	Control - no PHMB
MRSA	3	11
P. aeruginosa	3	10
E. cloacae	0	2
S. aureus	0	1

Reference: Product Review: Using PHMB antimicrobial to prevent wound infection by Keith Moore, David Gray - Wounds UK, 2007, Vol 3, No 2