

Introducing

Silicone Sacrum Foam Dressing with Adhesive Border

Safe n Simple focuses on product quality, patient safety, and cost containment. Simpurity wound care dressings follow evidence based medicine, facilitates moist wound healing and improves patients outcomes while lowering costs.



FEATURES:

- Comfortable and conformable, easy to apply and remove, has repositionable adhesive
- High absorbent foam pad with gentle silicone adhesive border
- Provides a moist environment for wound healing
- Minimizes trauma to the wound and reduces pain during dressing changes

Product Details:

Size	Item#	Box	Case	HCPC
7"x7"	SNS77777	12 Box	16 Boxes	A6213

INDICATIONS:

Simpurity Foam Wound Dressings are indicated for the management of moderately-to-highly exuding wounds such as:

- Pressure sores of the sacral area
- Chronic and acute exudative sacral wounds
- Lacerations and abrasions on shallow wounds

CONTRAINDICATIONS:

Simpurity Foam Wound Dressings are not indicated for use on the following:

- Third-degree burns
- Dry wounds
- Surgical implantations

DRESSING STRUCTURE:

Simpurity Silicone Foam Wound Dressings consist of a multi-layer design construction for maximum protection

